
Cod.: P303MED050

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1. Display LED
2. Vano batterie
3. Sensore infrarossi
4. Cordino

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## ENGLISH - FINGERTIP PULSE OXIMETER

### Description of product Fig.A

1. LED Display
2. Battery case
3. Infrared sensor
4. Hanging cord

### General description

Oxygen saturation is a blood index that reflects the percentage of hemoglobin saturated with oxygen compared to the total amount of hemoglobin present in the blood, and it is a very important physiological parameter in relation to respiratory and circulatory systems. A reduction of oxygen saturation in blood can be caused by various respiratory diseases, but also by imbalances linked to particular clinical conditions (anesthesia, trauma, invasive clinical exams, etc.). In these conditions, it is thus very important to monitor the level of oxygen saturation in the blood.

### Diagram of operation Fig.B

1. Infrared-ray receiving tube
2. Infrared-ray transmitting tube

### Scope of application

This oximeter is a portable, non-invasive device that can be used to measure the oxygen saturation of arterial hemoglobin (SpO2) and heart rate both in adult and children patients. This device cannot be used for continuous monitoring. This device does not require any routine calibration or other maintenance besides replacing of batteries.

### General safety rules

Do not use the oximeter when the patient is undergoing MRI or CT exams. Explosion hazard: do not use the oximeter in explosive environments. The oximeter is only intended as an adjunct in patient clinical assessment. It must be used in conjunction with other clinical evaluation methods of symptoms and conditions. Frequently check the position of oximeter sensor to determine the correct positioning in relation to the circulation and sensitivity of the patient's skin. Read this manual carefully before using the device. The oximeter is not equipped with SpO2 alarms, so it not designed for continuous monitoring.

Change the position of the sensor and verify the integrity of the skin, the circulation and the correct alignment every 2 hours at least.

Inaccurate measurements may be caused by autoclaving sanitizing process, or ethylene oxide sterilizing, or by immersing the device in liquids. Significant levels of dysfunctional hemoglobin (such as carboxyl-hemoglobin or methemoglobin) may affect the measurement.

Intravascular dyes such as indocyanine green or methylene blue can cause inaccurate readings.

SpO2 measurement can be negatively affected in the presence of high ambient light. Cover the sensor (for instance with a surgical towel) if necessary.

Sudden movements of patient can cause inaccurate readings.

Interferences caused by high frequency electrosurgical units or interferences causes by defibrillators can cause inaccurate readings.

Venous pulsations can cause unprecise measurements.

Do not place the sensor on a limb where a blood pressure monitor cuff, an arterial catheter or intravascular line is placed.

Do not use the device in case of hypotension, severe vasoconstriction, severe anemia or hypothermia.

Do not use the device in case of cardiac arrest or shock. Bright nails, painted nails or lacquered nails can cause inaccurate readings or missed readings.

Do not immerse the device in water. Do not use the device in environments where flammable anesthetic mixtures with air, oxygen or nitrous oxide are kept.

### Display Fig.C

1. Oxygen saturation
2. Pulse rate
3. SpO2 wave shape

### Installing of batteries

Insert 2 type AAA batteries inside the battery case, following the correct polarity. Close the battery cover by sliding it horizontally following the arrow.

**Important: insert the batteries following the correct polarity, to avoid damaging the device.**

Remove batteries if the oximeter is not used for a long period of time.

### Installing of hanging cord Fig.D

Insert the narrowest part of the cord through the eyelet and place the end part until the end of the cord before pulling.

## Operation and measurement

After correctly installing the batteries, open the clamp and insert the finger with the nail facing upwards. Fig. E

After a few seconds, display will turn on, and the device will produce a beeping sound. The display will show the signal related to the wave shape of the saturation level, and after a few seconds the values related to the oxygen saturation and pulse rate will appear.

**Important: keep the finger steady during measuring; it is recommended to remain seated or to lay down during measuring.**

Once measuring is done, remove the finger from the device, that will automatically shut off when signal is absent.

**Important: when display shows a low battery level, replace batteries to ensure correct operation of the device and accuracy of measurement.**

After each use, clean the rubber part inside the oximeter that comes in contact with finger with alcohol (rubber inside the device is medical rubber, it does not contain any toxin, and it is safe for the skin).

## Technical specifications

Display type: OLED

SpO2:

measuring range 70%-100%

precision 80%~99%, ±2% (included 80%); 70%~79%, ±3%; < 70% no requirement.

Pulse rate:

measuring range 30-254BPM.

Parameters of LED probe:

	Wave length	Radiation power
RED	660 ± 2nm	1.8mW
IR	940 ± 10nm	2.0mW

Power: two alkaline batteries type AAA.

Auto stand by mode: device will automatically shut off when it does not detect any signal (in approx. 8 seconds)

Size: 58x36x33 mm

Operation environmental requirements:

Temperature 5°~40°C

Humidity: 15%-80%

Atmospheric pressure: 70kPa~106kPa

Transport and storage environmental requirements:

Temperature -10°~40°C

Humidity: 10%-95%

Atmospheric pressure: 50kPa~106kPa

Non-corrosive gas and well ventilated environment

## Cleaning and maintenance

Replace batteries when the indicator on the display shows low battery level.

Before and after use, clean the rubber part inside the oximeter that comes in contact with finger with alcohol

Remove batteries from battery case if the device is not used for a long period of time.

Keep the device in a place with room temperature between -10°~40° (14~104°F) and humidity 10%-95%.

This device must be always stored in a dry place. A humid environment can damage the device reducing its life time.

**Warning: do not immerse the device in water or other liquids.**

## Troubleshooting

Problem	Possible cause	Solution
SpO2 or PR values cannot be shown normally	1. Finger in not placed deep enough inside the oximeter. 2. Patient's oxyhemoglobin level is too low to be detected	1. Place the finger correctly and try again. 2. Try several times. If you are sure the oximeter is working correctly, please seek medical assistance immediately.
SpO2 or PR values detected are unsteady	1. Finger in not placed deep enough inside the oximeter. 2. Finger is not steady or the patient moved	1. Place the finger correctly and try again. 2. Try again remaining still

Problem	Possible cause	Solution
The device does not turn on	1. Batteries are dead or not inserted. 2. Batteries are not inserted correctly. 3. The oximeter is damaged.	1. Replace batteries. 2. Insert batteries again. 3. Check with your local distributor.
Display indicators turn off suddenly	1. Device is programmed to automatically shut off when no signal is detected for more than 8 seconds. 2. Batteries are low.	1. Normal situation 2. Replace batteries.

## Symbols and definitions

	Impermeability index
	Manufacturer
	Lotto
	Read this manual carefully before use
	BF type device
	EU representative
	CE certificate

**For any improvement reasons, Beper reserves the right to modify or improve the product without any notice.**

 The European directive 2011/65/EU on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream.

Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment. The crossed out "wheeled bin" symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected.

## GUARANTEE CERTIFICATE

This appliance has been checked in the factory. From the date of original purchase a 24 month guarantee applies to material and production defects. The purchase receipt and the guarantee certificate must be submitted together in case of claim to guarantee. **The guarantee is valid only with a guaranty certificate and proof of purchase (fiscal receipt) indicating the date of purchase and the model of the appliance**

For any technical assistance, please contact directly the seller or our head office in order to preserve the appliance's efficiency and NOT TO void the guarantee. Any intervention on this appliance by non-authorized persons will automatically void the guarantee. Contact your distributor in your country or after sales department beper.

E-mail [assistenza@beper.com](mailto:assistenza@beper.com) which will forward your enquiries to your distributor.

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#### Troubleshooting

Problem	Possible cause	Solution
SpO2 or PR values cannot be shown normally	1. Finger in not placed deep enough inside the oximeter. 2. Patient's oxyhemoglobin level is too low to be detected	1. Place the finger correctly and try again. 2. Try several times. If you are sure the oximeter level is working correctly, please seek medical assistance immediately.
SpO2 or PR values detected are unsteady	1. Finger in not placed deep enough inside the oximeter. 2. Finger is not steady or patient moved	1. Place the finger correctly and try again. 2. Try again remaining still
The device does not turn on	1. Batteries are dead or not inserted. 2. Batteries are not inserted correctly. 3. The oximeter is damaged.	1. Replace batteries. 2. Insert batteries again. 3. Check with your local distributor.
Display indicators turn off suddenly	1. Device is programmed to automatically shut off when no signal is detected for more than 8 seconds. 2. Batteries are low.	1. Normal situation 2. Replace batteries.

#### Symbols and definitions

IP22	Impermeability index
	Manufacturer
	Lotto
	Read this manual carefully before use
	BF type device
	EU representative
	CE certificate

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© 2023
**LOT B022** Distrib: Beper S.r.l. - beper.com
Xuzhou Yongkang Electronic Science Technology Co.,Ltd
Profnrx GmbH