Cod.: P303MED100

MISURATORE DI PRESSIONE DA BRACCIO	IT
ARM BLOOD PRESSURE MONITOR	EN
ARMDRUCKMESSGERÄT	DE
TENSIOMÈTRE BRAS	FR
TENSIÓMETRO DE BRAZO	ES
TLAKOMJER ZA NADLAKTICU	HR
MERILNIK KRVNEGA TLAKA ZA NADLAKET	SL
ΠΙΕΣΟΜΕΤΡΟ ΒΡΑΧΙΟΝΑ	GR
TENSIOMETRU PENTRU BRAȚ	RO
TLAKOMĚR NA PAŽI	CZ
ROKAS ASINSSPIEDIENA MĒRĪTĀJS	LV
ŽASTINIS KRAUJOSPŪDŽIO MATUOKLIS	LT
ÖLAVARRE VERERÖHUMÖÖTJA	EE
BLOEDDRUKMETER VOOR DE BOVENARM	NL
АПАРАТ ЗА КРЪВНО НАЛЯГАНЕ ЗА РЪКА	BG
MERAČ KRVNOG PRITISKA ZA NADLAKTICU	RS
FELKAROS VÉRNYOMÁSMÉRŐ	HU

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Safety notice

For specific information about your blood pressure, please consult your doctor.

Explanation of symbols:

	Please read this manual thoroughly before using the unit. Please retain this manual for future reference. Operate the unit only as intended.
	Caution: Do not mix new and old batteries simultaneously. Do not mix batteries of different brands. Install batteries paying attention to polarity indicated. Remove batteries from the device in case you will not use for long time. Replace batteries when on the display will appear this symbol  Dispose batteries properly, observe local laws and regulation.
	Appliance type B
	Factory
	Forbidden: Do not allow children and / or persons with reduced capacity to intend to use this unit. Do not attempt to repair the product yourself. When using this product do not use mobile phones or other devices that may generate electromagnetic fields.
	Mandatory: Contact your doctor if the results are regularly abnormal. Use the appliance only as described.
	Authorized representative in the European community
	This electrical device must not be disposed of with household waste.
	The product complies with the requirements of Regulation (EU) 2017/745 (MDR) on medical devices.
	Date of manufacture
	Medical device
	Unique device identifier
IP21	Device protected against the ingress of solid objects with a diameter of 12 mm or greater, and against condensation droplets. Suitable for indoor use.

Product Description Fig. A

1. Air tube connector
2. WHO (World Health Organization) blood pressure classification indicator
3. Backlit LCD display
4. Systolic pressure indicator
5. Diastolic pressure indicator
6. Pulse rate indicator
7. Mode button
8. Power button
9. Set button
10. Battery compartment
11. Cuff
12. Air tube connector
13. Air tube
14. USB Type-C cable port

Display description Fig. B

1. Groups
2. Low battery indicator
3. WHO (World Health Organization) blood pressure classification indicator
4. Pressure detection indicators
5. Average of last 3 readings
6. Irregular heartbeat indicator
7. Time
8. Date
9. Pulse rate
10. Diastolic pressure
11. Systolic pressure
12. Heartbeat indicator
13. Number of stored tests

Some important guidelines

Read this manual carefully before using the device and keep it for future reference. For specific information regarding blood pressure, consult your physician.

The digital blood pressure monitor is intended for measuring blood pressure and pulse rate on the upper arm. Use this product only for its intended purpose. Use of this equipment outside the specified environment may damage the equipment and result in inaccurate measurements.

This product is intended for use by adults and adolescents over 12 years of age.

This appliance must not be used by persons (including children) with reduced physical, sensory, or mental capabilities, or by persons lacking experience and knowledge of the appliance, unless they are supervised or instructed concerning the use of the appliance by a person responsible for their safety.

Avoid using this equipment adjacent to or stacked with other equipment, as this may result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Arm Blood Pressure Monitor Use Instructions

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) to any part of the fully automatic digital blood pressure monitor. Failure to maintain this separation distance may result in degradation of the equipment's performance. In the event of loss or degradation of essential performance caused by electromagnetic disturbances, the operator should stop using the system and contact customer service.

The use of accessories or cables other than those specified may result in increased electromagnetic emissions or decreased electromagnetic immunity, which could impair the proper operation of the equipment.

The user of this product is the patient and is considered the intended operator. The patient must not perform any maintenance operations other than battery replacement. Do not modify this equipment without authorization. If calibration or operational deviations are suspected, contact the service center. The heat resistance of the device will be maintained throughout the expected service life of the medical electrical equipment.

Do not confuse blood pressure monitoring with self-diagnosis. Measurement values should only be interpreted by a physician or healthcare professional familiar with the user's medical history.

Consult your physician if repeated measurements indicate abnormal or unusual values.

Individuals with severe circulatory disorders may experience discomfort during measurement. Consult your physician before use.

This blood pressure monitor is not recommended for individuals with severe arrhythmias.

People suffering from specific circulatory disorders should use the blood pressure monitor only after consulting their physician. In such cases, the oscillometric measurement method may provide inaccurate readings.

Excessively frequent measurements may cause injury or discomfort to the patient due to temporary interference with blood circulation. Do not apply the cuff over wounds or skin lesions, as this may aggravate the condition.

Fully unroll the connection tube before use. If compressed or bent, the connection tube may cause continuous pressure in the cuff, resulting in interference with blood flow and possible injury to the patient.

Use only the cuff approved for this device. The use of other cuffs may result in incorrect measurement readings. Do not use the cuff if it is damaged.

Keep the device out of reach of children and infants. This product contains small parts that may present a choking hazard if swallowed. Do not leave the cuff unattended. It may cause suffocation if swallowed by a child or wrapped around the neck.

The cuff and blood pressure monitor housing have been tested for biocompatibility and do not contain allergenic or harmful materials. Discontinue use if an allergic reaction occurs during operation.

Avoid eating, exercising, or bathing within 30 minutes before taking a measurement.

Sit in a quiet environment for at least 5 minutes before taking a measurement.

Do not stand during measurement.

During measurement, sit in a relaxed position and keep your arm at heart level.

Avoid talking or moving any part of your body during measurement.

Wait at least 3 minutes before taking another measurement.

It is recommended to measure blood pressure at the same time of day, on the same arm, and preferably in the same position.

Battery warnings

This product is designed to operate with the specified type of batteries, not included. Use only the type and polarity specified in the manual or on the device. Do not attempt to recharge non-rechargeable batteries, and do not disassemble, crush, short-circuit, or dispose of batteries in fire. Remove batteries if the product will not be used for an extended period and dispose of them according to local regulations. Damaged or leaking batteries can cause burns or irritation if they come into contact with skin or eyes; avoid direct contact with leaked liquid. Always keep batteries out of reach of children, as small batteries can pose a choking hazard or cause serious internal injury if swallowed. If you suspect ingestion or introduction of batteries into any part of the body, seek immediate medical attention.

Power supply

The product operates either with four AAA alkaline batteries (not included) or connected to a USB source using a USB Type-C cable (not included). If using batteries, slide the battery compartment cover as indicated by the arrow. Install four new AAA alkaline batteries, observing the correct polarity. Close the battery compartment cover.

If operating via USB, before first use, remove the protective film from the USB Type C port on the side of the product and insert the USB cable, then insert the USB plug into a USB power socket. This blood pressure monitor works with standard 5V – 1000mA (1A) USB adapters available on the market. You can use a medical-grade USB adapter with double insulation (2MOPP), which meets the highest standards for medical devices.

System settings

With the product turned off but powered by batteries or connected to a USB port, press the “SET” button to access the system settings. The display lights up and the memory group icon starts flashing. In the system settings, test results can be stored in two different groups. This allows multiple users to save individual test results (up to 75 memories per group). Press the “MODE” button to select the group. Test results will be automatically stored in each selected group.

Once the group is selected, press the “SET” button again to set the time/date mode. First set the hour by pressing the “MODE” button. Press the “SET” button again to confirm the hour and then move on to the minutes. Continue setting the month and day in the same way. Each press of the “SET” button saves the setting and proceeds in sequence (hour, minutes, month, day).

After finally selecting the day, press “SET” again to display the memory group at the top; at this point, you can press “MODE” to move on to the next group and proceed with the settings in the same way.

Settings saved

In any setting mode, press the “ON/OFF” button to turn off the unit.

All information will be saved.

If the unit is left on and unused for three minutes, all information will be automatically saved and the device will turn off.

Use

After installing the batteries or connecting the device via USB cable, and with the system set up, apply the cuff by inserting the connector on the plastic tube of the cuff into the appropriate slot on the main unit. (Fig. C)

Remove any clothing from the arm.

Arm Blood Pressure Monitor Use Instructions

Rest for several minutes before testing. Sit in a quiet place, preferably at a desk or table, with your arm resting on a stable surface and your feet on the floor. (Fig. D)

Apply the cuff to your left arm and hold it at heart level. Wrap the cuff around your arm without over-tightening it (you can use about a finger's width between the cuff and your arm as a guide). The bottom of the cuff should be positioned approximately 1-2 cm above the elbow. Fig. E

At this point, you can proceed with the test.

Press the "ON/OFF" button until you hear a beep. The display lights up and a quick diagnosis appears for one second. A long beep indicates that the device is ready for measurement.

 **Warning: The unit will not function if residual air from the previous test is present in the cuff.**

The blood pressure monitor begins inflating and measuring blood pressure. Remain relaxed during the test. Avoid talking or moving other parts of your body.

During the measurement, the display remains illuminated and an arrow flashes until the blood pressure stabilizes.

The initial pressure is first pumped to 190 mmHg. If the user's current systolic blood pressure is above 190 mmHg, the unit automatically reinflates to the appropriate level.

 **Warning: Pressurization will gradually decrease and eventually stop if the cuff is not properly applied to the arm. In this case, press the "ON/OFF" button to turn off the unit.**

After the cuff is inflated, the air slowly decreases as indicated by the cuff pressure reading. A heart symbol flashes simultaneously on the display to indicate a heartbeat.


 **Warning: If the cuff pressure becomes too extreme during testing, press the "ON/OFF" button to immediately turn off the device. The cuff pressure will quickly dissipate once the unit is turned off.**

Four short beeps will sound when the test is complete.

The screen displays systolic and diastolic blood pressure readings.

An indicator representing the current reading is displayed next to the corresponding WHO classification.

The product features a classification indicator based on guidelines established by the World Health Organization. The graph below indicates the test results.

After reading the results, turn off the product by pressing the ON/OFF button. Otherwise, if no button is pressed, the product will automatically turn off after 3 minutes.

Warning: Make sure you have selected the appropriate memory group before testing.

If the number of tests exceeds the 75 assigned memories per group, the most recent tests will appear first, thus eliminating the oldest readings.

Warning: Do not be alarmed by an abnormal reading. An individual's blood pressure is determined by 2-3 readings taken at the same time each day over an extended period of time. Consult your doctor if test results remain abnormal.

If the symbol  appears on the display and four beeps sound, the batteries are low and should be replaced.

If the word "Err" appears on the display, the cuff may not be positioned correctly or the blood pressure may be above 300 mmHg. Try adjusting the cuff and/or relaxing, then try again.

Irregular heartbeat indicator

If the monitor detects an irregular heart rhythm two or more times during the measurement process, a heart symbol, representing an irregular heartbeat, will appear alongside the measurement results. Irregular heart rhythm is defined as a heart rate that is 25% slower or faster than the average heart rate detected during systolic and diastolic blood pressure measurements. Consult your doctor if the irregular heartbeat symbol appears frequently during the test results.

Average of the last 3 tests

With the device turned off, press the "MODE" button to activate the on-screen display. The unit performs a self-test; the screen will display the average of the last 3 readings for the last group used. The "AVG" symbol will appear along with the corresponding WHO blood pressure indicator. The memory check mode can be accessed by pressing the "MODE" buttons. To check the average results for other groups, select the desired group before pressing the "MODE" button with the device turned off.

Viewing and clearing memory

With the device turned off, you can check past test results using the “MODE” button. When you press the button, the device emits a beep and the display lights up. The most recent and oldest test results in memory can be viewed by repeatedly pressing and holding the “MODE” button. You can scroll through all stored test results. The test numbers are displayed in succession at the bottom right of the display, followed by the blood pressure readings.

 **Warning:** Previous test results will only be displayed from the most recently used memory group. To check previous test results from other memory groups, you must first select the desired group and then turn off the monitor.

The memory of the selected group can be cleared in memory check mode. Press and hold the “SET” button for about 3 seconds to delete all memory results from the selected group. The monitor emits four short beeps to indicate deletion and then goes directly into test mode. Press the “ON/OFF” button to turn off the unit.

 **Warning:** Once cleared, memory cannot be recovered.

WHO Measurement Classification Indicator

The blood pressure monitor is equipped with a blood pressure classification indicator based on the guidelines established by the World Health Organization (WHO).

WHO classification	Systolic (mmHg)	Diastolic (mmHg)
Hypertension grade 3	≥180	≥110
Hypertension grade 2	160-179	100-109
Hypertension grade 1	140-159	90-99
High – normal	130-139	85-89
Normal	120-129	80-84
Optimal	<120	<80

Maintenance

Avoid dropping, knocking, or throwing the unit.

Avoid extreme temperatures. Do not expose the unit to direct sunlight.

The equipment can be cleaned by a lay operator following the procedures described below.

Before cleaning the device, ensure that it is disconnected from the USB power supply.

To clean the unit, use a soft cloth slightly dampened with water.

To clean the cuff, apply a small amount of alcohol to a soft cloth. Avoid allowing any liquid to enter the blood pressure monitor. Dry all parts thoroughly before storing the product.

 **Warning:** Do not immerse the cuff in water.

Remove the batteries if not in use for an extended period of time.
It is recommended to check the performance at least every 2 years.

Specifications

Product description	Upper Arm Blood Pressure Monitor
Model number	P303MED100 (DBP-13C2)
Display	Digital LCD - Size: 66mm x 43mm
Measuring method	Oscillometric
Measurement intervals	Blood Pressure: 0 mmHg - 300 mmHg Pulse Rate: 30 to 180 beats per minute
Measurement accuracy	Blood Pressure: ± 3 mmHg Pulse Rate: $\pm 5\%$
Pressurization	Automatic
Memory	2x150 memories with date and time
Functions	Low battery indicator WHO classification indicator Auto power off Irregular heartbeat indicator Average of last 3 readings Backlight
Power supply	4 AAA batteries or USB Type-C cable DC 5.0V 1A
Batteries life	Approx. 2 months based on 3 measurements per day
Unit weight	Approx. 310g (excluding batteries and USB cable)
Dimensions	Approx. 139 x 88 x 43mm
Cuff Circumference	Adjustable circumference 22 - 42 cm
Protection Rating	IP21
Use Environment	Temperature: 10°C - 40°C (50°F - 104°F) Humidity: 15% - 90% RH Atmospheric pressure: 700hPa - 1060hPa
Storage Environment	Temperature: -25°C - 70°C Humidity: $\leq 93\%$ RH

Directives and standards

This product meets the requirements of the EU standard for non-invasive medical devices. It has been certified according to EC directives and bears the CE mark (mark of conformity) "CE 0123." The product complies with European standards EN 60601-1, EN 60601-1-2, EN 60601-1-11, EN 60601-1-6, ISO 81060-2, EN 80601-2-30, EN 62366-1, and EN 62304. The device meets the requirements of Directive (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

Electric compatibility information

Table 1

Guidance and declaration of manufacturer-electromagnetic emissions	
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.	
Emissions test	Compliance
Radiated emission CISPR 11	Group 1, class B.
Conducted emission CISPR 11	Group 1, class B.
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies

Table 2

Guidance and declaration of manufacturer-electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
IMMUNITY test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrostatic transient/burst IEC 61000-4-4	± 2 kV, 100 kHz, for AC power port	± 2 kV, 100 kHz, for AC power port
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV (differential mode)	± 0.5 kV, ± 1 kV (differential mode)
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50 Hz or 60 Hz	30 A/m; 50 Hz and 60 Hz
Radiated RF EM fields IEC 61000-4-3	3 V/m or 10 V/m 80 MHz-2.7 GHz 80% AM at 1kHz	10 V/m 80MHz-2.7 GHz 80%AM at 1kHz
Conducted disturbances Induced by RF fields IEC 61000-4-6	3 V in 0.15 MHz-80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz	3 V in 0.15 MHz-80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz

IMMUNITY test	IEC 60601 test level	Compliance level
Proximity magnetic fields IEC 61000-4-39	See Table 4	See Table 4

Table 3

Guidance and declaration of manufacturer-electromagnetic immunity				
Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. Arm-type Fully Automatic Digital Blood Pressure Monitor has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014+AMD1:2020. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and this medical equipment and/or systems as recommended below.				
Test frequency (MHz)	Band a) (MHz)	Service a)	Modulation	IMMUNITY TEST LEVEL (V/M)
385	380 a 390	TETRA 400	Pulse modulation b) 18 Hz	27
450	430 a 470	GMRS 460, FRS 460	FM c) ±5 kHz deviation 1 kHz sine	28
710 745 780	704 a 787	LTE Band 13,17	Pulse modulation b) 217 Hz	9
810 870 930	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	28
1720 1845 1970	1700 to 1990	GSM 1800, CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS	Pulse modulation b) 217 Hz	28

2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	28
5240 5500 5785	5100 to 5800	WLAN 802.11a/n	Pulse modulation b) 217 Hz	9
If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3. a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50 % duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case.				

Table 4

Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields		
Test frequency	Modulation	IMMUNITY TEST LEVEL (A/m)
30 kHz a)	CW	8
134,2 kHz	Pulse Modulation b) 2,1 kHz	65 c)
13,56 MHz	Pulse Modulation b) 50 kHz	7,5 c)
a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT. b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) r.m.s., before modulation is applied.		

In the pursuit of continuous improvement, Beper reserves the right to make modifications and improvements to the product without prior notice.

The European directive 2011/65/EU on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment. The crossed out "wheeled bin" symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected.

GUARANTEE CERTIFICATE

This appliance has been checked in the factory. From the date of original purchase a 24 month guarantee applies to material and production defects. The purchase receipt and the guarantee certificate must be submitted together in case of claim to guarantee.

The guarantee is valid only with a guaranty certificate and proof of purchase (fiscal receipt) indicating the date of purchase and the model of the appliance

For any technical assistance, please contact directly the seller or our head office in order to preserve the appliance's efficiency and NOT TO void the guarantee. Any intervention on this appliance by non-authorized persons will automatically void the guarantee.

GUARANTEE CONDITIONS

If the appliance shows defects as a result of faulty material and/or production during the period of guarantee, we guarantee repair free of charge on condition that:

- The appliance has been used properly and for the purpose for which it has been intended.
- The appliance has not been tampered with, otherwise it is not maintainable.
- The purchase receipt shall be presented.
- The appliance showing a fair wear and tear shall not be covered by this guarantee.

Therefore, any part that could be accidentally broken or having visible signs of use in consumable products (such as lamps, batteries, heating elements...) the aesthetic parts are excluded from the guarantee, and any defect whatsoever resulting from non-respect of the rules for use, negligence in usage and/or maintenance of the appliance, carelessness, wrong or improper installation, damage during transport and any other damage not attributable the supplier.

For every defect that could not be repaired within the guarantee period, the appliance will be replaced free of charge.

In any case, if the part to be replaced for defect, breakage or malfunction is an accessory and/or a detachable part of the product, Beper reserves the right to replace only the very part in question and not the entire product

Contact your distributor in your country or after sales department Beper.

E-mail assistenza@beper.com which will forward your enquiries to your distributor.

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*Sign up for the newsletter to receive
the latest Beper news in advance*

Digital Ear Thermometer

**SCAN
ME!**

beper.com 

CUSTOMER CARE BEPER

BEPER SRL

Via Salieri, 30 - 37050 - Vallese di Oppeano - Verona
Tel. 045/7134674 – Fax 045/6984019
beper.com