Cod.: 40.500

MASSAGGIATORE ANTICELLULITE	IT
ANTI-CELLULITE MASSAGER	EN
MASSEUR ANTI-CELLULITE	FR
ANTI-CELLULITE-MASSAGEGERÄT	DE
MASAJEADOR ANTICELULÍTICO	ES
ΣΥΣΚΕΥΗ ΜΑΣΑΖ ΚΑΤΑ ΤΗΣ ΚΥΤΤΑΡΙΤΙΔΑΣ	GR
APARAT DE MASAJ ANTICELULITIC	RO
LIPOMASÁŽ NA CELULITIDU	CZ
ANTI-CELLULITE MASSAGE APPARAAT	NL
PROFESJONALNY SYSTEM LIPOMODELINGOWY	PL
ANTI-CELULĪTA MASAŽIERIS	LV
TSELLULIIDIIVASTANE MASSEERIJA	EE
ANTICELULIT MASAŽER	RS
ANTICELULITNI MASAŽNI APARAT	SL
ANTICELIULITINIS MASAŽUOKLIS	LT

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General warnings

Read the following instructions before using the appliance.

Before and during use of the appliance, some basic precautions must be followed.

After removing all packaging materials, check the integrity of the appliance. In case of any doubt do not use the appliance and contact professionally qualified personnel. Packaging materials (plastic bags, styrofoam, etc.) must always be kept out of the reach of children, as potential cause of risk.

Always make sure that the mains voltage equals to the voltage indicated on the technical data label and that the electrical system is compatible with the power of the appliance.

Never unplug the appliance from the socket by pulling the power cord. Make sure the cable never gets in contact with hot or sharp surfaces. Do not use the appliance if the power cord is damaged.

If the power cord is damaged, it must be replaced by the manufacturer or by its Technical After Sales service or by a qualified person, in order to avoid any potential risk.

Only connect the appliance to an AC power outlet.

It is generally not advisable to use adapters, multiple sockets and/or cord extensions.

When their use is necessary only use adapter and cord extensions compliant with the applicable safety regulations.

This appliance must only be used for the operation for which it was expressly designed. Any other use is to be considered incorrect and thus dangerous, causing the expiration of warranty. The manufacturer cannot be held responsible for any damage caused by incorrect, improper and unreasonable use.

To avoid dangerous overheating completely unwind the power cord and unplug the socket from the power outlet when the appliance is not in use. Before carrying out any cleaning or maintenance operation, unplug the socket from the power outlet.

Do not immerse the appliance in water or other liquids.

Keep the appliance away from heating sources (e.g. radiator).

This appliance must not be used by persons (including children) with limited physical, sensory or mental capabilities, by persons that lack experience and knowledge of the appliance, unless they are closely watched or well instructed by a person responsible for their safety regarding the use of the appliance.

Children must not play with the appliance.

When the appliance must be disposed of, it is recommended to make it inoperative.

It is also recommended to render harmless those parts of the appliance likely to constitute a danger. This appliance must not be used by children. Children must not play with the appliance. Keep the appliance and its cable away from children under 8 years of age.

This appliance can be used by children of 8 years or older and by persons with limited physical, sensory or mental capabilities, or by persons that lack experience and knowledge of the appliance, provided they have been given complete instructions concerning the safe use of the appliance and provided that they understand the related risks.

Warnings

Read this manual carefully, it contains important instructions for a safe installation, use and maintenance of the appliance.

This unit can be used by children aged 8 and above and individuals with reduced physical, sensory or mental capacity, or lack of experience or knowledge, if they are provided with adequate supervision or have received instructions regarding the use of safety equipment and that they understand the risks related.

Cleaning and maintenance should not be carried out by children without supervision.

If the power cable gets damaged, it has to be replaced only at a service center approved by the manufacturer.

In the presence of medical conditions, such as phlebitis, varicose veins, skin excoriations or other issues, consult your doctor for frequency and duration of applications, as well as possible contraindications.

Expectant mothers, people who suffer from diabetes and people who carry pacemakers must consult their doctor before using this appliance.

Individuals insensible to heat must be careful when using the appliance.

Do not use the massager for more than 15 minutes.

Before making a new cycle of massages, allow the appliance to cool down.

Do not massage the same body area for more than 5 minutes.

After removing the packaging, check the appliance for any signs of damage. In case of visible damage, do not use it and contact a qualified technician.

If the plug and socket are not compatible, the socket has to be replaced with a suitable type by a qualified technician.

When not in use, disconnect the appliance from the power supply and make sure that it is switched off. Never touch the appliance with wet or humid hands.

Do not leave the appliance exposed to weathering (rain, sunlight)

Do not drop or hit the appliance.

In case of a breakdown or malfunctioning of the appliance, disconnect it at once. Do not attempt to repair it, but take it to a service center approved by the manufacturer.

The appliance is not designed for medical purposes.

Exercise extreme care when using the massager on your neck and back area. Always make sure to undo your hair before massaging the neck and back area, to prevent the movement of the heads from accidentally trapping some strands of hair, causing them to tear.

Description Fig.A

1. Main unit
2. Power/massage/massage + IR switch
3. Adjustable massage intensity controller
4. Infrared light
5. Massage head adaptor
6. Cord
7. Finger massage head
8. Spiky roller massage head
9. Ball massage head
10. Wavy surface massage head
11. Protective cloth head

Tips for use

Self-massage is best done in the morning because the body is rested and more receptive. During the evening, especially after a long day spent sitting down or on your feet, legs are already congested and additional stress may increase the inflammatory state, so it is best to apply only cream with a light massage at this time of day.

Massage the areas to be treated in small concentric circles. The friction helps to improve local circulation. Move slowly and fluidly, starting from the bottom up and without applying excessive pressure to the tissues.

Massaging the legs and feet is especially good for alleviating the sensation of tiredness and heaviness and promoting a pleasant relaxed feeling.

Use instructions

Choose the desired accessory for the massage, mount it onto the hexagonal base making sure that the appliance is both switched off and unplugged. Fig.B

Plug the appliance and place the switch to 1 position, to activate the massage only, or 2 to also activate the infrared light. Adjust the intensity of the massage with the controller.

The accessories supplied can be used interchangeably.

Check the diagram in Fig. C and D to choose the speed and the right accessory for every body area to be treated.

Use the cloth head if you are applying cream during the massage to prevent dirtying the accessories. Mount it by turning it clockwise onto the base unit, making sure that it is properly closed.

Accessories 7 and 9 are known as multinodes due to the various shaped protrusions. Using these, you can perform a percussion massage of varying intensity. Applying greater pressure will result in a deeper massage.

The roller accessory 8 is used to perform a vibration massage that works against cellulite thanks to the up and down movement. It also promotes good blood circulation.

Accessory 10 is particularly ideal for a more delicate and less intense friction massage. It is recommended for alleviating muscular pain from cramps due to the mild, relaxing action of the massage.

The infrared/heat function combined with vibration massage applies heat to the skin which helps absorb the nutrients from the cosmetic products used.

Do not massage the same area for more than 5 minutes and do not run the appliance for more than 15 minutes at a time. Before beginning a new session of massages wait at least 15 minutes until the appliance cools down.

Cleaning and maintenance

Unplug the appliance before cleaning it.

Do not use solvents, detergents or abrasive products.

Use a dry cloth to clean the outer surface.

When the appliance is not in use, unplug it.

To prevent the cable from becoming damaged, do not twist it, pull it or wind it around the appliance. If it is damaged, have it replaced immediately by qualified personnel.

 **Warning: do not immerse the appliance in water or other liquids.**

Technical data

Power: 25W

Power supply: 230V~ 50Hz

For any improvement reasons, Beper reserves the right to modify or improve the product without any notice.

The European directive 2011/65/EU on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment.

The crossed out "wheeled bin" symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected.

GUARANTEE CERTIFICATE

This appliance has been checked in the factory. From the date of original purchase a 24 month guarantee applies to material and production defects. The purchase receipt and the guarantee certificate must be submitted together in case of claim to guarantee.

The guarantee is valid only with a guaranty certificate and proof of purchase (fiscal receipt) indicating the date of purchase and the model of the appliance

For any technical assistance, please contact directly the seller or our head office in order to preserve the appliance's efficiency and NOT TO void the guarantee. Any intervention on this appliance by non-authorized persons will automatically void the guarantee.

GUARANTEE CONDITIONS

If the appliance shows defects as a result of faulty material and/or production during the period of guarantee, we guarantee repair free of charge on condition that:

- The appliance has been used properly and for the purpose for which it has been intended.
- The appliance has not been tampered with, otherwise it is not maintainable.
- The purchase receipt shall be presented.
- The appliance showing a fair wear and tear shall not be covered by this guarantee.

Therefore, any part that could be accidentally broken or having visible signs of use in consumable products (such as lamps, batteries, heating elements...) the aesthetic parts are excluded from the guarantee, and any defect whatsoever resulting from non-respect of the rules for use, negligence in usage and/or maintenance of the appliance, carelessness, wrong or improper installation, damage during transport and any other damage not attributable the supplier.

For every defect that could not be repaired within the guarantee period, the appliance will be replaced free of charge.

In any case, if the part to be replaced for defect, breakage or malfunction is an accessory and/or a detachable part of the product, Beper reserves the right to replace only the very part in question and not the entire product

Contact your distributor in your country or after sales department Beper.

E-mail assistenza@beper.com which will forward your enquiries to your distributor.

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Neck massager

Cod.: P302MAS100

EMS technology (low frequency impulse)
Flexible electrodes perfectly adaptable to the neck
EMS massage with adjustable intensity
6 EMS massaging modes and heating function
Lithium battery
USB charging cable included

A woman with long brown hair is smiling and wearing the teal and grey neck massager. The device is shown in use, highlighting its fit and the control dial.

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USB

YOU MIGHT ALSO LIKE

Rechargeable Percussion Massager

Cod.: P302MAS002

- 6 adjustable intensity levels
- 4 interchangeable massaging accessories
- Long lasting lithium battery
- Ergonomic and anti-slip handle
- Low noise and light weight for easy use

BEPER SRL

Via Salieri, 30 - 37050 - Vallese di Oppeano - Verona
Tel. 045/7134674 – Fax 045/6984019

beper.com

CUSTOMER CARE BEPER