# *beper*

IT	BILANCIA PESAPERSONE IMPEDENZIOMETRICA	pag. 2
EN	IMPEDANCE BATHROOM SCALE	pag. 5
FR	PÈSE-PERSONNE IMPÉDANCEMÈTRE	pag. 8
DE	IMPEDANZ-PERSONENWAAGE	pag. 11
ES	BÁSCULA DE BIOIMPEDANCIA	pag. 14
GR	ΖΥΓΑΡΙΑ ΜΠΑΝΙΟΥ ΒΙΟΗΛΕΚΤΡΙΚΗΣ ΑΝΤΙΣΤΑΣΗΣ	pag. 17
RO	CÂNTAR DE DIAGNOSTIC	pag. 20
CZ	IMPEDANČNÍ OSOBNÍ VÁHA	pag. 23
NL	IMPEDANTIE BADKAMER WEEGSCHAAL	pag. 26
PL	WAGA DIAGNOSTYCZNA	pag. 29
LV	VANNAS ISTABAS SVARI	pag. 32

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**assistenza@beper.com**

## **General safety advice**

### **Read the operating instructions carefully before first use of the appliance.**

Read the following instructions carefully, since they give you useful safety information about installation, use and maintenance and help you avoid mishaps and possible accidents.

Remove the packing and make sure that the appliance is intact, with special attention to the supply cable. The elements of the packing (plastic bags, polystyrene, etc) do not have to be accessible to the children, in order to preserve them from danger sources.

It is recommended not to throw these packing material in the domestic waste, but to deliver them to the appropriate station of collection or destruction refusals, asking eventually information to the manager of the sanitation department of your own municipality.

This appliance is to be destined to its special purpose only.

Any other use is to be considered as improper and thus dangerous. The manufacturer cannot be responsible for any damages deriving from improper, wrong or incautious use.

To avoid any heating risk, it is recommended to unroll completely the supply cord and to unplug the appliance, when not used.

Before cleaning, always unplug the appliance.

To protect against electric shock, do not immerse cord, plug or any parts of the appliance into water or other fluids.

Do not hold the appliance with wet hands or feet.

Do not hold the appliance with naked feet.

Do not expose the appliance to atmospheric agents such as rain, wind, snow, hail.

Do not use the appliance near explosives, highly flammable materials, gases, burning flames, heaters.

Do not allow children, disabled people or every people without experience or technical ability to use the appliance if they are not properly watched over. Close supervision is necessary when any appliance used by or near children.

Make sure that children do not play with this appliance.

If you decide not to use the appliance any longer, make it unserviceable by simply cutting the supply cable after disconnecting the plug from the socket.

**Also make any dangerous parts unserviceable, especially for children who could play with the appliance.**

## Some warnings

In case of a fault and/or malfunction, turn off and do not tamper with the scale. For any repairs, contact only a centre authorized by the manufacturer and request the use of original spare parts. Failure to comply with the above may compromise security and invalidate the warranty.

This scale is for domestic use only, not for use in hospital diagnosis. The use of the scale is not indicated for people with implanted medical equipment, such as pacemakers.

Do not leave the scale exposed to atmospheric agents (rain, sun, etc.).

Do not use the scale too much close to radio or TV to avoid electromagnetic interference.

Do not immerse the scale in water.

Do not place the scale too much close to heat sources.

Do not place the scale in very humid places.

Do not place the scale on a wet surface.

Only place the scale on a flat and stable surface, do not place it on a carpet as it would not be stable.

Do not leave the scale upright.

Do not place objects on the scale.

Do not drop the scale or subject it to severe shocks.

Do not expose the scale to excessive weights.

Do not jump on the scale.

Do not allow liquids to enter the inside of the scale or spray water or cleaning detergents directly on the scale.

Before using the scale, remove the transparent film in the battery compartment, if present.

Use only the type of batteries specified in the instruction manual.

**⚠ Important:** when measuring body mass and water, make sure you are barefoot, clean and completely in contact with the plant on the surface of the scale to avoid percentage errors in body mass.

### **⚠ Caution**

Discharged batteries become dangerous for the possible loss of electrolyte fluid that is very corrosive. To avoid damage to the appliance caused by this possibility, replace the exhausted batteries without delay. Remember to remove the batteries when you are not going to use the balance for a long time.

Damages caused by battery liquid are not recognized.

### **Battery replacement**

When the scale does not turn on or displays "Lo" it means that the batteries are exhausted. Replace them with new ones.

Carefully, turn the balance upside down, so that the bottom of the scale is accessible.

Press the battery compartment in the direction indicated by the arrow and open the cover by tapping it.

Insert the batteries into the compartment according to the correct polarity. Incorrect installation can damage the scale. Close the cover and then turn the scale back on.

### Operation

Remove the packing and place the scale on a solid and flat surface.

Before using the scale, remove the plastic sheet from the battery compartment to enable it.

Press "SET" button and wait till "0.0" will be shown on the display.

Get on the balance with both feet, hold them in the glass surface as much as possible and stay still for a few seconds.

The display will automatically display the corresponding weight automatically.

If the display shows "Err" it means there was a weighing error and then we need to do it again. Get off the balance and repeat.

The scale automatically shuts off after weighing.

Weighing up to 180 kg. To avoid any damage, if the "EEEE" message appears on the display immediately get off from the scale as it indicates it is beyond the weighing limit.

Body weight is subject to changes during the day, so there may be weight differences. Lightweight weight changes are usually normal. It is advisable to always weigh yourself at the same time, preferably before breakfast.

The balance is a precise tool and is very sensitive to movement of the body. For better precision: climb on the balance, relax your body, hold your head firm and stay still.

**Note: there is a small switch near the battery compartment to switch the weight from kg to lb or st.**

### Set personal data

The scale has the optional function of weighing and finding out your water and body fat percentage. To do this you must first input the following data in the scale memory: gender, age and height. The scale memory can store data up to 10 people.

To input the data in the memory, follow the next steps:

Press and release the SET button and wait for the display shows " 0.0".

Press the SET button and use UP/DOWN button to find the memory number under which you wish to save your personal data (10 memories available: P01-P10). Then press the SET button to confirm your choice.

Select male or female by using UP/DOWN buttons and confirming by pressing SET button.

Log in your age using the UP/DOWN buttons and confirming by using SET button.

Log in your height using the UP/DOWN buttons and confirming by using SET button. Note: you can also use the unit selector button under the scale to select height measure (Inches / Centimetre).

Once you have set your profile, press the SET button and select your profile. Your saved data will briefly flash onto the screen of the scale.

These steps can be repeated to program the remaining 9 memories if desired.

### Getting weighed and obtaining body fat, body water, bone and muscle percentage Fig.A

Once all the data have been introduced into the memory, you can obtain the calculation of your body mass.

Follow these steps:

Press and release the SET button and wait for the digits "0.0" appear on the display. Press SET again and use the UP " or DOWN button

until your profile number appears on the screen. Wait for the digits "0.0" appear on the display.

Place down the scale and place your body at the centre on top of the metal sensors. Remember to be bare feet and to stay still.

The scale will show flashing "000" while busy calculating. Once done, the scale will flash results in the following sequence:

1. Body Fat %
2. Body Water % (Hydration)
3. Bone %
4. Muscle %

The results will repeatedly flash for 3 times. After stepping down the scale will turn off automatically.

### Calculate body fat percentage

This appliance has the built-in option to calculate the percentage of your body fat by using bioelectrical impedance techniques. During calculation are also involved the following variables: height, age, gender, weight.

Optimum results according to standards generally applied to adults:

Men: 14-24%

Women: 20-30%

Any results below these indicators are considered underweight and higher than these are considered overweight. In the event of doubt, consult your doctor for clarification.

<b>Body fat ranges</b>							
<b>Male</b>	Age	10-20	20-30	30-40	40-50	50-60	60-99
	Underweight	< 10%	< 11%	< 13%	< 14%	< 15%	< 16%
	Healthy	10-18%	11-19%	13-21%	14-22%	15-23%	16-24%
	Overweight	>18%	>19%	>21%	>22%	>23%	>24%
<b>Female</b>	Age	10-20	20-30	30-40	40-50	50-60	60-99
	Underweight	< 15%	< 15%	< 17%	< 18%	< 19%	< 19%
	Healthy	15-23%	15-23,5%	17-25%	18-26%	19-27%	19-27%

The scale will also use one of the following descriptions as per the results:

Underweight

Healthy

Overweight

Obese

#### **Cleaning the scale**

The scale can be cleaned with a slightly dampened soft cloth with warm water.

Never clean with diluents, abrasive products in general; these could damage the surface of the scale.

Do not immerse the scale in water.

Protect from liquid infiltration.

Keep away the scale from preparation for hair, cosmetics or bath solutions; they can cause damage to the aesthetic parts and internal scale circuit.

#### **Technical features**

Maximum capacity: 180Kg

Power supply: lithium 1xCR2032 3V (included)

**For any improvement reasons, Beper reserves the right to modify or improve the product without any notice.**

The European directive 2011/65/EU on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment.

The crossed out "wheeled bin" symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected.

### **GUARANTEE CERTIFICATE**

This appliance has been checked in the factory. From the date of original purchase a 24 month guarantee applies to material and production defects. The purchase receipt and the guarantee certificate must be submitted together in case of claim to guarantee.

**The guarantee is valid only with a guaranty certificate and proof of purchase (fiscal receipt) indicating the date of purchase and the model of the appliance**

For any technical assistance, please contact directly the seller or our head office in order to preserve the appliance's efficiency and NOT TO void the guarantee. Any intervention on this appliance by non-authorized persons will automatically void the guarantee.

### **GUARANTEE CONDITIONS**

If the appliance shows defects as a result of faulty material and/or production during the period of guarantee, we guarantee repair free of charge on condition that:

- The appliance has been used properly and for the purpose for which it has been intended.
- The appliance has not been tampered with, otherwise it is not maintainable.
- The purchase receipt shall be presented.
- The appliance showing a fair wear and tear shall not be covered by this guarantee.

Therefore, any part that could be accidentally broken or having visible signs of use in consumable products (such as lamps, batteries, heating elements...) the aesthetic parts are excluded from the guarantee, and any defect whatsoever resulting from non-respect of the rules for use, negligence in usage and/or maintenance of the appliance, carelessness, wrong or improper installation, damage during transport and any other damage not attributable the supplier.

For every defect that could not be repaired within the guarantee period, the appliance will be replaced free of charge.

In any case, if the part to be replaced for defect, breakage or malfunction is an accessory and/or a detachable part of the product, Beper reserves the right to replace only the very part in question and not the entire product

Contact your distributor in your country or after sales department Beper.

E-mail [assistenza@beper.com](mailto:assistenza@beper.com) which will forward your enquiries to your distributor.

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