beper

- MASSAGGIATORE A INFRAROSSI - MANUALE DI ISTRUZIONI
- BODY INFRARED MASSAGER - USE INSTRUCTIONS
- MASSEUR Á INFRAROUGES - MANUEL D'INSTRUCTIONS
- INFRAROT MASSAGE-GERÄT - BETRIEBSANLEITUNG
- MASAJEADOR CON INFRARROJOS - MANUAL DE INSTRUCCIONES
- ΣΥΣΚΕΥΗ ΜΑΣΑΖ ΜΕ ΥΠΕΡΥΘΡΕΣ - ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ

Cod.: 40.942

ITALIANO	pag. 4
ENGLISH	pag. 9
FRANÇAIS	pag. 14
DEUTSCH	pag. 20
ESPAÑOL	pag. 26
ΕΛΛΗΝΙΚΗ	pag. 32

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e-mail: assistenza@beper.com

GENERAL SAFETY ADVICE

READ THE OPERATING INSTRUCTIONS CAREFULLY BEFORE FIRST USE OF THE APPLIANCE.

Read the following instructions carefully, since they give you useful safety information about installation, use and maintenance and help you avoid mishaps and possible accidents.

Remove the packing and make sure that the appliance is intact, with special attention to the supply cable. The elements of the packing (plastic bags, polystyrene, etc) do not have to be accessible to the children, in order to preserve them from danger sources.

It is recommended not to throw these packing material in the domestic waste, but to deliver them to the appropriate station of collection or destruction refusals, asking eventually information to the manager of the Sanitation Department of your own municipality.

Before connecting the appliance make sure that the mains voltage corresponds to the data specified on the label.

If the plug does not suit the socket, let it be replaced with another suitable plug by qualified personnel only, which shall also ensure that the section of the socket cables is suitable for your appliance power absorption.

Never pull the supply cord to unplug the appliance.

Make sure that the supply cord never comes into contact with hot or cutting surfaces

Do not use this appliance if the supply cord is damaged. In case of damage, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.

Only connect this appliance on alternating current

Generally, the use adaptors multiple plugs and/or extensions is not recommended; should they be necessary, please use only simple or multiple adaptors and extensions in compliance with the safety rules in force, making sure that the capacity absorption limit marked on the simple adaptors and the extensions and the maximum power limit marked on the multiple adaptor are not exceeded.

This appliance is to be destined to its special purpose only.

Any other use is to be considered as improper and thus dangerous. The manufacturer cannot be responsible for any damages deriving from improper, wrong or incautious use.

To avoid any heating risk, it is recommended to unroll completely the supply cord and to unplug the appliance, when not used.

Before cleaning, always unplug the appliance.

To protect against electric shock, do not immerse cord, plug or any parts of the appliance into water or other fluids.

Do not hold the appliance with wet hands or feet.

Do not hold the appliance with naked feet.

Do not expose the appliance to atmospheric agents such as rain, wind, snow, hail.

Do not use the appliance near explosives, highly flammable materials, gases, burning flames, heaters.

Do not allow children, disabled people or every people without experience or technical ability to use the appliance if they are not properly watched over. Close supervision is necessary when any appliance used by or near children.

Make sure that children do not play with this appliance.

If you decide not to use the appliance any longer, make it unserviceable by simply cutting the supply cable after disconnecting the plug from the socket.

ALSO MAKE ANY DANGEROUS PARTS UNSERVICEABLE, ESPECIALLY FOR CHILDREN WHO COULD PLAY WITH THE APPLIANCE.

USE ADVICE

Make sure that the mains voltage corresponds to the data specified on the label of the cord supply. The massage cushion is to be destined to its special purpose, only for household use and cannot be used for commercial purpose.

Always unplug the appliance from the socket when changing accessories or cleaning.

Do not let the infrared massage without supervision when plugged to the socket.

Keep the cord supply away from hot surfaces.

Keep away from children, except if there is a close supervision by a responsible person.

Do not use the infrared massage near to eyes or every delicate parts of body.

Regularly check the massage cushion, do not use it if damaged.

Do not try to open the infrared massage or parts of it.

Do not use the infrared massage in damp places; keep away from heat sources or sun. Do not use outdoor.

Do not immerge the infrared massage into water or other liquids. If the infrared massage is wet, unplug immediately it from the socket.

Do not use the infrared massage with damp or wet hands.

Do not use the infrared massage under a blanket or in a close room ; if these rules are not respected there is a risk of fire or injuries.

Do not use the infrared massage if damaged; contact the manufacturer or authorized centre.

If the infrared massage is not used, unplug the cord from the socket.

Do not let the infrared massage without supervision if there are children. This appliance cannot be used by children. Keep away from children.

Do not cover the appliance when working.

SECURITY PRECAUTION

Use the infrared massage only for its purpose. This appliance cannot be used for medical reason and cannot substitute medical cure. Only for pleasant massage.

If any doubt on your health, contact your doctor before using the massage cushion.

People wearing a pacemaker or pregnant women must be visited by a doctor before using this appliance.

The infrared massage must not be used by people with physical problems that could prevent them from a correct use of the control panels.

Do not use the infrared massage on children, disabled people or asleep ; do not use if you are affected by loss of the sensitivity or blood circulation problems without being visited by a doctor.

Do never use the infrared massage directly if ignition or cutaneous eruption on the skin. Do not use the infrared massage before going to sleep, the massage has a stimulant effect that could not make you sleeping.

When using the cushion massage if you feel unpleasant feelings, stop the appliance and contact your doctor. Do not use for longer periods than those recommended, use it according to the particular physical characteristics of each person. This infrared massage is not a professional appliance and has been manufactured only for personal use, in order to practise a beneficial massage on tired muscles. Do not use the infrared massage instead of medical care. Do not use the infrared massage for more than 30 minutes. The excessive use can cause breakdown. If that happens, stop immediately the appliance and wait until the appliance is completely cool to use it again.

PRODUCT DESCRIPTION Fig. 1

1. Housing of interchangeable heads
2. Head massage
3. Head specialized massage (deep)
4. Exfoliating head
5. Infrared function LED indicators
6. Body massager
7. Massage intensity adjustment
8. Operating mode selector
(off - vibrating massage action -
action vibrating massage with infrared)
9. Ergonomic handle
10. Power Cord

INSTRUCTION Fig.2

The massage cushion is perfect to eliminate muscular or physical tiredness, the pains caused by a bad blood circulation or muscular tension. The massage cushion has been studied for the massage of each part of the body, a comfortable way and with effects similar to those of a manual massage. Through the action of vibrating head you get an immediate improvement in blood circulation and muscle tissue, relief from fatigue and reduce the production of lactic acid.

Select the desired head, insert it into the slot :

Massage head to eliminate muscle fatigue, tiredness or just relaxing

Massage head for a deep massage to loosen muscular tension

Exfoliating head for a smooth and soft skin

Plug the power cord (10) into the electrical outlet.

Select the desired action using the appropriate operating mode selector. **Fig.3**

L = vibrating massage action with infrared

0 = off

I = vibrating massage action

The function infrared can accelerate blood circulation, not hot.

The massager is more effective if you use it without clothes, it's a personal and individual choice, can also be used with clothing.

Use the massager for a varying period from person to person according to their physical characteristics, if necessary, consult a doctor.

CLEANING THE MASSAGER

After each use, always switch off the massager by placing the selector switch to 0 and unplug the power cord from the electrical socket.

Clean the body (6) and heads (2-3-4) with a soft dry cloth.

Do not immerse in water or other liquids.

TECHNICAL DATA

Power 25 W

Power supply: 230V ~ 50Hz

For any improvement reasons, Beper reserves the right to modify or improve the product without any notice.

The European directive 2011/65/EU on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment. The crossed out "wheeled bin" symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected.

GUARANTEE CERTIFICATE

This appliance has been checked in the factory. From the date of original purchase a 24 month guarantee applies to material and production defects. The purchase receipt and the guarantee certificate must be submitted together in case of claim to guarantee.

The guarantee is valid only with a guaranty certificate and proof of purchase (fiscal receipt) indicating the date of purchase and the model of the appliance

For any technical assistance, please contact directly the seller or our head office in order to preserve the appliance's efficiency and NOT TO void the guarantee. Any intervention on this appliance by non-authorized persons will automatically void the guarantee.

GUARANTEE CONDITIONS

If the appliance shows defects as a result of faulty material and/or production during the period of guarantee, we guarantee repair free of charge on condition that:

- The appliance has been used properly and for the purpose for which it has been intended.
- The appliance has not been tampered with, otherwise it is not maintainable.
- The purchase receipt shall be presented.
- The appliance showing a fair wear and tear shall not be covered by this guarantee.

Therefore, any part that could be accidentally broken or having visible signs of use in consumable products (such as lamps, batteries, heating elements...) the aesthetic parts are excluded from the guarantee, and any defect whatsoever resulting from non-respect of the rules for use, negligence in usage and/or maintenance of the appliance, carelessness, wrong or improper installation, damage during transport and any other damage not attributable the supplier.

For every defect that could not be repaired within the guarantee period, the appliance will be replaced free of charge.

In any case, if the part to be replaced for defect, breakage or malfunction is an accessory and/or a detachable part of the product, Beper reserves the right to replace only the very part in question and not the entire product

TECHNICAL ASSISTANCE

Even after the guarantee period we will always pay attention to repairs for defected appliances. For technical assistance and/or repairs after the guarantee period you can directly contact the below address:

CONTACT YOUR DISTRIBUTOR IN YOUR COUNTRY OR ATER SALES DEPARTMENT BEPER.
E-MAIL assistenza@beper.com WHICH WILL FORWARD YOUR ENQUIRIES TO YOUR DISTRIBUTOR.

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