beper

IT	MASSAGGIATORE PIEDI MULTIFUNZIONE	pag. 2
EN	MULTIFUNCTIONAL FOOT MASSAGER	pag. 8
FR	MASSEUR DE PIEDS MULTIFONCTIONNEL	pag. 14
DE	MULTIFUNKTIONALES FUSSMASSAGEGERÄT	pag. 21
ES	MASAJEADOR DE PIES MULTIFUNCIÓN	pag. 27
GR	ΠΟΛΥΛΕΙΤΟΥΡΓΙΚΗ ΣΥΣΚΕΥΗ ΜΑΣΑΖ ΠΟΔΙΩΝ	pag. 33
RO	APARAT DE MASAJ PENTRU PICIOARE MULTIFUNCȚIONAL	pag. 40
CZ	MULTIFUNKČNÍ MASÁŽNÍ PŘÍSTROJ NA NOHY	pag. 46
NL	MULTIFUNCTIONEEL VOETMASSAGEAPPARAAT	pag. 52
LV	DAUDZFUNKCIONĀLS PĒDU MASIERIS	pag. 58
SL	VEČNAMENSKI MASAŽER ZA STOPALA	pag. 64

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assistenza@beper.com

Precautions

Read the following instructions carefully, since they give you useful safety information about installation, use and maintenance and help you avoid mishaps and possible accidents.

Remove the packing and make sure that the appliance is intact, with special attention to the supply cable.

The elements of the packing (plastic bags, polystyrene, etc) do not have to be accessible to the children, in order to preserve them from danger sources. It is recommended not to throw these packing material in the domestic waste, but to deliver them to the appropriate station of collection or destruction refusals, asking eventually information to the manager of the Sanitation Department of your own municipality.

If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.

Before connecting the appliance make sure that the mains voltage corresponds to the data specified on the label.

The appliance can be damaged if the voltage is not correct.

If the plug does not suit the socket, let it be replaced with another suitable plug by qualified personnel only, which shall also ensure that the section of the socket cables is suitable for your appliance power absorption.

Generally, the use adaptors multiple plugs and/or extensions is not recommended; should they be necessary, please use only simple or multiple adaptors and extensions in compliance with the safety rules in force, making sure that the capacity absorption limit marked on the simple adaptors and the extensions and the maximum power limit marked on the multiple adaptor are not exceeded.

Do not use the appliance if the cable or the plug shows any damages.

Do not use the appliance after a possible irregular working. In this case, switch it off and do not break it open. Refer to an authorized Beper after-sale service for repair and require original spare parts only. Failing to follow these instructions may impair the safety of your appliance.

This appliance is to be destined to its special purpose only. Any other use is to be considered as improper and thus dangerous. The manufacturer cannot be responsible for any damages deriving from improper, wrong or incautious use.

This appliance is for household use only.

To protect against electric shock, do not immerse cord, plug or any parts of the appliance into water or other fluids.

Do not use the appliance near explosives, highly flammable materials, gases or burning flames.

Do not allow children, disabled people or every people without experience or technical ability to use the appliance if they are not properly watched over. Close supervision is necessary when any appliance used by or near children.

Warnings for use

Read this manual carefully before use, and keep it for future reference.

Do not leave the appliance unattended when connected to the wall socket; unplug it when not in use, and before removing the foot sleeves.

Remove the adapter from the appliance if not in use for long periods of time.

Do not use the appliance under blankets or pillows, to avoid over heating and consequent risks of fire, electric shock or injuries.

Do not use this appliance for self-treatment of particular medical conditions that require the attention of specialized medical personnel.

Do not use the appliance in case of malaise.

This massager is to be considered a non-professional appliance, designed to relieve muscular pain and soreness. It is not to be considered as a substitute for specialized medical care. If symptoms and ailments should not disappear or should worsen, please consult your doctor.

⚠ Warning: consult with specialized medical personnel is required before use of the appliance in case of pregnancy or post-partum, cardiovascular and circulatory problems, pacemaker, diabetes or other particular conditions, phlebitis and/or thrombosis, joint dysfunction, inflammation or injuries, weak bones, osteoporosis, high fever, scoliosis or other similar conditions, prosthesis, or other medical conditions and diseases that require specialized medical attention.

Do not use the appliance while sleeping or feeling sleepy.

Do not use the appliance in case of severely swollen, or inflamed feet, or in case of skin rash and inflammation.

The appliance should not be used by individuals with reduced skin sensitivity or with circulatory problems.

Immediately stop using the appliance in case of pain or discomfort related to the use of the massager.

The appliance produces therapeutic heat: individuals with high sensitivity to heat should exercise extreme care during use of the appliance.

Do not use this appliance in connection with other therapeutic devices, electric blankets, etc.

Use this appliance only for its intended use, and as described in this manual. Do not use any accessory/component not provided with the appliance.

Do not use the appliance under the influence of alcoholic beverages, pain killers or sedatives.

Do not let the children play with the appliance.

This device cannot be used by children, or persons with reduced physical, sensitive or intellectual capabilities, or with lack of experience and knowledge, without supervision or if they are not given instructions and warnings on the correct and safe use of the device and the potential risks related to its use.

Do not transport the appliance by holding it from the power cable and do not use the cable as a handle.

Do not leave the appliance near heating sources, such as radiators, stoves, or other devices that generate heat.

Make sure that the vent holes are always free from dust and dirt.

Do not drop the appliance and do not insert any object inside the appliance.

Do not operate the appliance in case the power cable or its plug are damaged, or in case of malfunctioning, accidental drop, or in case the appliance has been accidentally dropped in water.

Do not use the appliance in the presence of aerosol products or where oxygen is being administered.

Do not stand while using the massager, only use it while sitting down.

Do not use the appliance if it has been immersed in water. Unplug it from the wall socket immediately.

Keep the power cable and the adapter away from heating sources.

It is advisable not to use the appliance one hour before and after meals, to avoid potential adverse effects.

The remote control provided contains a pre-installed lithium battery. Keep the remote control and its battery away from children.

⚠ Warning: if swallowed, the battery can cause internal burns, perforation of soft tissues, and lead to potentially lethal injuries. Severe burns can occur within 2 hours from swallowing. Seek immediate medical attention if ingestion is suspected.

Do not let children replace the battery without adult supervision.

Do not expose battery to high temperature. Remove it from remote control if it is not used for long time.

Do not dispose of the exhausted battery in regular waste; do not throw it in fire or attempt to burn it. Follow applicable local regulations for the correct disposal of batteries.

Description of product Fig. A

1. Upper cover
2. Control panel
3. Zipper for removal of foot sleeves
4. Detachable foot sleeves
5. Lower cover
6. Adapter
7. Remote control

Operation from control panel Fig. B

Connect the adapter to the cable found at the base of the appliance, making sure the connectors fit perfectly, and insert the plug into the wall socket: the LED indicator of the POWER key starts flickering, indicating that the appliance is in stand-by mode.

Pressing the POWER key once, the massager starts operating in WORK mode (P1): low heat, low intensity.

At this point, press the MODE key to select the desired massage program:

P1 – shiatsu massage + air compression – green LED light

P2 – air compression massage – blue LED light

P3 – shiatsu massage – red LED light (in this mode, a light compression will be activated, for a better massaging experience).

Once the desired program has been selected, it is possible to set the level of air compression and heat.

Press the HEAT key to adjust the heat on low level (green LED light) or high level (red LED light). Press the INTENSIFY key to adjust air compression on 3 different intensity levels: low intensity (green LED light), medium intensity (blue LED light), high intensity (red LED light).

At the beginning, it is better to set the air compression on low intensity, and then increase it gradually according to your needs (when moving from one level to another, it takes about 2 minutes to fully feel the difference in intensity).

⚠ Warning: especially at the beginning of the massage, it is normal to feel a little pain and, at times, numbness, due to the intensity of the massage and the work on the foot muscles and joints; if the pain or the discomfort persist and/or it should become unbearable, stop using the appliance immediately.

⚠ Warning: the over heating protection system fitted in the appliance limits the heating speed of the product; it can take about 5-10 minutes before the heat can be perceived noticeably.

Operation from remote control Fig.B

Activate the battery by removing the plastic protection strip.

All the functions available from the control panel can also be activated from the remote control; refer to previous paragraph.

Auto shut off

For safety reasons, this massager is equipped with an auto shut off system, which stops the operation of the appliance after 15 minutes of use, and bringing it to stand by mode (LED indicator light on POWER key flickers). Press POWER key to resume operation of the massager, and select the desired function.

Cleaning and maintenance

Before cleaning the appliance, make sure to disconnect the adapter from the wall socket.

Use the zipper to remove the foot sleeves.

The foot sleeves can be washed in washing machine at 30°C, mild fine wash (do not tumble dry and do not iron). Use a soft cloth, slightly dampened with water, to clean the outer surface; if necessary, use a neutral detergent to remove stains. Dry thoroughly.

⚠ Warning: do not use benzene or other abrasive or corrosive detergents.

When not in use, store the massager away from dust, humidity and direct sunlight, preferably inside its original packaging.

Technical data

Power: 36W

Adapter:

Input: AC 100-240V, 50/60Hz, 0.8A

Output: DC 24V, 1.5A

For any improvement reasons, Inside reserves the right to modify or improve the product without any notice.

The European directive 2011/65/EU on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment. The crossed out “wheeled bin” symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected.

GUARANTEE CERTIFICATE

This appliance has been checked in the factory. From the date of original purchase a 24 month guarantee applies to material and production defects. The purchase receipt and the guarantee certificate must be submitted together in case of claim to guarantee.

The guarantee is valid only with a guaranty certificate and proof of purchase (fiscal receipt) indicating the date of purchase and the model of the appliance

For any technical assistance, please contact directly the seller or our head office in order to preserve the appliance's efficiency and NOT TO void the guarantee. Any intervention on this appliance by non-authorized persons will automatically void the guarantee.

GUARANTEE CONDITIONS

If the appliance shows defects as a result of faulty material and/or production during the period of guarantee, we guarantee repair free of charge on condition that:

- The appliance has been used properly and for the purpose for which it has been intended.
- The appliance has not been tampered with, otherwise it is not maintainable.
- The purchase receipt shall be presented.
- The appliance showing a fair wear and tear shall not be covered by this guarantee.

Therefore, any part that could be accidentally broken or having visible signs of use in consumable products (such as lamps, batteries, heating elements...) the aesthetic parts are excluded from the guarantee, and any defect whatsoever resulting from non-respect of the rules for use, negligence in usage and/or maintenance of the appliance, carelessness, wrong or improper installation, damage during transport and any other damage not attributable the supplier.

For every defect that could not be repaired within the guarantee period, the appliance will be replaced free of charge.

In any case, if the part to be replaced for defect, breakage or malfunction is an accessory and/or a detachable part of the product, Beper reserves the right to replace only the very part in question and not the entire product

Contact your distributor in your country or after sales department beper.

E-mail assistenza@beper.com which will forward your enquiries to your distributor.

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