ITALIANO

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e-mail: assistenza@beper.com - **www.beper.com**

ENGLISH
Warning :
Read instructions thoroughly before using digital thermometer
Choking Hazard: Thermometer cap and battery may be fatal if swallowed. Do not allow children to use this device without parental supervision.
Do not use thermometer in ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
Do not place thermometer battery near extreme heat as it may explode.
Note: Use of the probe cover may result in a 0.1°C (0.2°F) discrepancy from actual temperature.
It is recommended the performance should be checked every two years.
Remove battery from the device when not in operation for a long time.
The performance of thermometer may be degraded if one or more of the following occur:
a) Operation outside stated temperature and humidity range;
b) Storage outside stated temperature and humidity range;
c) Mechanical shock;
d) Patient temperature is below ambient temperature.
This digital thermometer provides a quick and highly accurate reading of an individual's body temperature. To better understand its functions and to provide years of dependable results, please read all instructions first.
This thermometer is manufactured according to the system of EN12470-3:2000+A1:2009, ASTM E1112-00.

PRECAUTION
* Do not bend or drop thermometer. It is not shockproof.
* Do not boil probe to disinfect.
* Do not store in direct sunlight or in areas of extreme heat.
* Do not disassemble. See BATTERY REPLACEMENT to replace battery. Do not put in trash.
*Recycle or manage as hazardous waste. Battery contains mercury.
*Avoid strong electromagnetic interference such as microwave ovens and cell phones.

→ Direct Current .
⤴ Type B Applied Part
⚠ Caution, consult accompanying documents
LOT Batch code.
🏢 Manufacturer
🇪🇺 European Authorized Representative
DIRECTIONS
1. Wipe or dip thermometer probe with isopropyl rubbing alcohol to disinfect .
2. Press the On/Off Button next to LCD display. A tone will sound as the screen shows ▯▯▯▯ followed by last recorded temperature.

After showing the self-test temperature, the thermometer is now in the testing mode.

***Note:** If the measured temperature is less than 32°C or 90°F , the LCD will display Lo, if the measured temperature is higher than 42.9°C or 109.9°F , the LCD will display Hi.

3. Position thermometer in desired location (mouth, rectum, or armpit.)

a) **Oral Use:** Place thermometer under tongue as indicated by position shown in Figure 2. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. Normal temperature between 35.7°C and 37.3°C (96.3°F and 99.1°F)

b) **Rectal Use:** Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1cm (less than 1/2) into rectum. Normal temperature between 36.2°C and 37.7°C (97.2°F and 99.9 °F).

c) **Armpit Use:** Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side.

From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required. Normal temperature between 35.2°C and 36.7°C (95.4°F and 98.1 °F).

4. The degree sign flashes throughout the testing process. When flashing stops an alarm will beep for approximately 10 seconds. The measured reading will appear on the LCD simultaneously. The displayed temperature will not change when thermometer is removed from its testing position

*Note: Normally Alarm beeps m ore rapidly w hen temperature reaches 37.8°C (100°F) or higher ,and the buzzes are" Bi- Bi-Bi- Bi-"; he buzzes are" Bi-Bi-Bi----- Bi-Bi-Bi----- Bi-Bi-Bi"

5. To prolong battery life, press the On/Off Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around10 minutes.

***C/ °F SWITCHABLE**

Temperature readings are available in the Fahrenheit or Celsius scale (°C/°F; located in the upper right corner of LCD.) With the unit off, press and hold the On/Off Button for approximately 4 seconds to change the current setting.

BATTERY REPLACEMENT

- Replace battery when appears in the lower right corner of LCD display.
- Pull battery cover off
- Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 1.5 DC button type LR41 or SR41,UCC392 , or equivalent. Be sure batte-ry is installed with polarity facing up.
- Slide battery chamber back into place and attach cover.

CLEANING AND DISINFECTION

Wipe the thermometer with a soft clean cloth.

For stubborn stains, wipe the thermometer with a cloth that has been dampened with water or a neutral detergent solution and then wring thoroughly. Finish by wiping with a soft dry cloth.

Observe the following to prevent damage to the thermometer.

- Do not use benzene, thinner, gasoline or other strong solvents to clean the thermometer.
- Do not soak the sensing section in alcohol for long periods of time or attempt to sterilize it using hot water (water at a temperature of 50°C (122°F) or higher).
- Do not use ultrasonic washing to clean the thermometer.

TECHNICAL DATA

Measurement range: (32.0+42.9)°C/(90+109.9)°F

Power supply: 1.5 DC battery LR41 type

For any improvement reasons, Beper reserves the right to modify or improve the product without any notice.

GUARANTEE CERTIFICATE

This appliance has been checked in the factory. From the date of original purchase a 24 month guarantee applies to material and production defects. The purchase receipt and the guarantee certificate must be submitted together in case of claim to guarantee.

CONTACT YOUR DISTRIBUTOR IN YOUR COUNTRY OR AFTER SALES DE-PARTMENT BEPER. E-MAIL **assistenza@beper.com** WHICH WILL FORWARD YOUR ENQUIRIES TO YOUR DISTRIBUTOR.

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