Cod.: P303BIP051

BILANCIA PESAPERSONE IMPEDENZIOMETRICA	**IT**
IMPEDANCE BATHROOM SCALE	**EN**
IMPEDANZ-PERSONENWAAGE	**DE**
PÈSE-PERSONNE IMPÉDANCEMÈTRE	**FR**
BÁSCULA DE IMPEDANCIA CORPORAL	**ES**
IMPEDANCIJSKA KUPAONSKA VAGA	**HR**
OSEBNA TEHTNICA Z IMPEDANČNO ANALIZO	**SL**
ΖΥΓΑΡΙΑΜΠΑΝΙΟΥ ΒΙΟΗΛΕΚΤΡΙΚΗΣ ΑΝΤΙΣΤΑΣΗΣ	**GR**
CÂNTAR CORPORAL CU ANALIZĂ PRIN IMPEDANȚĂ	**RO**
IMPEDANČNÍ OSOBNÍ VÁHA	**CZ**
ĶERMEŅA SASTĀVA SVARI AR BIOIMPEDANCES ANALĪZI	**LV**
KŪNO SUDĖTIES SVARSTYKLĖS	**LT**
KEHAKOOSTISE ANALÜÜSIGA VANNITOAKAAL	**EE**
IMPEDANTIE BADKAMERWEEGSCHAAL	**NL**
ПЕРСОНАЛЕН КАНТАР С ИМПЕДАНСЕН АНАЛИЗ	**BG**
DIGITALNA VAGA SA ANALIZOM TELESNE MASE	**RS**
TESTÖSSZETÉTEL-ELEMZŐ SZEMÉLYMÉRLEG	**HU**

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General warnings

Read the following instructions before using the appliance.

Before and during use of the appliance, some basic precautions must be followed.

After removing all packaging materials, check the integrity of the appliance. In case of any doubt do not use the appliance and contact professionally qualified personnel.

This appliance must only be used for the operation for which it was expressly designed. Any other use is to be considered incorrect and thus dangerous, causing the expiration of warranty. The manufacturer cannot be held responsible for any damage caused by incorrect, improper and unreasonable use.

Do not leave the appliance exposed to atmospheric agents (rain, sun, etc.). Do not keep the appliance near heat sources (e.g. radiator).

This appliance must not be used by persons (including children) with limited physical, sensory or mental capabilities, by persons that lack experience and knowledge of the appliance, unless they are closely watched or well instructed by a person responsible for their safety regarding the use of the appliance.

THIS APPLIANCE MUST NOT BE USED BY CHILDREN.

Some warnings

In case of a fault and/or malfunction, turn off and do not tamper with the scale. For any repairs, contact only a centre authorized by the manufacturer and request the use of original spare parts. Failure to comply with the above may compromise security and invalidate the warranty.

This scale is for domestic use only, not for use in hospital diagnosis.

The use of the scale is not indicated for people with implanted medical equipment, such as pacemakers.

Do not leave the scale exposed to atmospheric agents (rain, sun, etc.).

Do not use the scale too much close to radio or TV to avoid electromagnetic interference.

Do not immerse the scale in water.

Do not place the scale too much close to heat sources.

Do not place the scale in very humid places.

Do not place the scale on a wet surface.
Only place the scale on a flat and stable surface, do not place it on a carpet as it would not be stable.
Do not leave the scale upright.
Do not place objects on the scale.
Do not drop the scale or subject it to severe shocks.
Do not expose the scale to excessive weights.
Do not jump on the scale.
Do not allow liquids to enter the inside of the scale or spray water or cleaning detergents directly on the scale.
Before using the scale, remove the transparent film in the battery compartment, if present.
Use only the type of batteries specified in the instruction manual.

Some warnings
In the event of a fault or malfunction, turn off the scale and do not tamper with it.
This scale is for home use only, not for hospital diagnostics.
The scale is not intended for use by people with implanted medical devices, such as pacemakers.
Do not leave the scale exposed to atmospheric agents (rain, sun, etc.).
Do not use the scale in the immediate vicinity of a radio or television to avoid electromagnetic interference.
Do not immerse the scale in water.
Do not place the scale in the immediate vicinity of heat sources.
Do not place the scale in very humid areas.
Do not place the scale on a wet surface.
Place the scale only on flat, stable surfaces; do not place it on carpets, as this would be unstable.
Do not leave the scale in an upright position.
Do not place objects on the scale.
Do not drop the scale or subject it to violent impacts.
Do not expose the scale to excessive weight.
Do not jump on the scale.
Do not allow liquids to enter the scale, nor spray water or cleaning agents directly onto the scale.
Before using the scale, remove the transparent plastic tab in the battery compartment, if present.

Use only the type of batteries specified in the instruction manual.

Product description Fig.A
1. Electrodes
2. Display
3. SET button
4. Parameter setting and configuration buttons
5. Measurement unit selection button
6. Battery compartment

Warnings for Use
Make sure to place the scale on a flat, stable surface.
Avoid exposing the scale to shock and/or violent impacts; do not expose it to excessive moisture.
Remove the batteries from the compartment if the scale will not be used for extended periods of time.
This scale is intended for home use only and not for use in a medical diagnostic setting. Consult your doctor if the readings are too high or too low. Use of the device is not recommended for people with electronic medical devices, such as pacemakers.

Installing and Replacing Batteries
This scale operates with 2 type AAA alkaline batteries, not included.
Insert the batteries inside the battery case, paying attention to the polarity indicated in the case.
The display always shows the battery status indicator: when the display shows “lo” indicator, replace the exhausted batteries with new ones.

Battery warnings
This product is designed to operate with the specified type of batteries, not included. Use only the type and polarity specified in the manual or on the device. Do not attempt to recharge non-rechargeable batteries, and do not disassemble, crush, short-circuit, or dispose of batteries in fire. Remove batteries if the product will not be used for an extended period and dispose of them according to local regulations. Damaged or leaking batteries can cause burns or irritation if they come into contact with skin or eyes; avoid direct contact with leaked liquid. Always keep batteries out of reach of children, as small batteries can pose a choking hazard or cause serious internal injury if swallowed. If you suspect ingestion or introduction of batteries into any part of the body, seek immediate medical attention.

Setting the scale Fig.B
Press the SET button on the platform: the number corresponding to the selected user will flash for a few seconds. Use the adjustment buttons to select the desired user number (from 0 to 9), then press the SET button to confirm. After confirming the user, the height in cm will appear: use the adjustment buttons to set the height (from 100 to 240 cm), and press the SET button to confirm. The Age display will then appear; use the adjustment buttons to set the selection (from 10 to 99), confirming with the SET button. Finally, the gender icon (male or female) will appear on the display; select the desired option using the adjustment buttons. Once all the personal user values have been set, the display will show the height and age data in sequence, then 0.0kg: the scale is now ready for use. Once all the users have been set, you can select the user simply by pressing the SET button with the scale off, and use the adjustment buttons to locate the desired user number: this will flash for a few seconds on the display, then the scale will show 0.0kg, and will be ready for use.

Using the scale Fig.B
If you simply want to measure your weight, simply step on the scale (barefoot): the scale turns on automatically, and the display shows your weight (the measurement ends when the weight on the display flashes for a few seconds).

⚠ **Warning: The scale’s maximum capacity is 180 kg; if you exceed this limit, the display shows “oL” and the scale turns off automatically.**

To measure your body mass index (BMI), before stepping on the scale, press the SET button and use the arrows to select your user number.

Once you have selected your user and the display shows 0.0 kg, step onto the platform slowly with bare feet, making sure your feet are clean and firmly planted on the steel electrodes. Remain as still as possible while the scale measures your weight and the values for body fat, bone mass, muscle mass, water percentage, calories, and BMI.
Once the measurement is complete, the display shows “ooo,” indicating that the scale is processing the data collected through the electrodes. After a few seconds, the display shows all the indicators related to the measured values:
BMI – Body Mass Index
CALO – Calories
MSL – Muscle Mass %
BON – Bone Mass %
FAT – Fat Mass %
WAT – Water %

⚠ **Warning: If there is a problem with data collection during the measurement, the display shows “Err” and the scale turns off automatically. Step off the scale and repeat the measurement, being careful to remain as still as possible and ensuring your feet are firmly planted on the electrodes.**

Cleaning and maintenance
It is advisable to clean the steel electrodes after each use to ensure proper functioning and maximum measurement accuracy.
Use a soft, slightly damp cloth to remove dust and any dirt from the glass platform and electrodes.

⚠ **Warning: Do not spray or pour abrasive or aggressive detergents directly onto the platform.**

Do not immerse the scale in water or other liquids; avoid allowing water or liquids to penetrate the base of the scale, as this could damage it.

Technical data
Power supply: 2xAAA alkaline batteries (not included)

In view of continuous improvement, Beper reserves the right to make changes and improvements to the product without prior notice.

The European directive 2011/65/EU on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment. The crossed out “wheeled bin” symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected.

GUARANTEE CERTIFICATE

This appliance has been checked in the factory. From the date of original purchase a 24 month guarantee applies to material and production defects. The purchase receipt and the guarantee certificate must be submitted together in case of claim to guarantee.

The guarantee is valid only with a guaranty certificate and proof of purchase (fiscal receipt) indicating the date of purchase and the model of the appliance

For any technical assistance, please contact directly the seller or our head office in order to preserve the appliance's efficiency and NOT TO void the guarantee. Any intervention on this appliance by non-authorized persons will automatically void the guarantee.

GUARANTEE CONDITIONS

If the appliance shows defects as a result of faulty material and/or production during the period of guarantee, we guarantee repair free of charge on condition that:

- The appliance has been used properly and for the purpose for which it has been intended.
- The appliance has not been tampered with, otherwise it is not maintainable.
- The purchase receipt shall be presented.
- The appliance showing a fair wear and tear shall not be covered by this guarantee.

Therefore, any part that could be accidentally broken or having visible signs of use in consumable products (such as lamps, batteries, heating elements…) the aesthetic parts are excluded from the guarantee, and any defect whatsoever resulting from non-respect of the rules for use, negligence in usage and/or maintenance of the appliance, carelessness, wrong or improper installation, damage during transport and any other damage not attributable the supplier.
For every defect that could not be repaired within the guarantee period, the appliance will be replaced free of charge. In any case, if the part to be replaced for defect, breakage or malfunction is an accessory and/or a detachable part of the product, Beper reserves the right to replace only the very part in question and not the entire product

Contact your distributor in your country or after sales department Beper.
E-mail assistenza@beper.com which will forward your enquiries to your distributor.

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